A-178

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EN Lumbar Belt A-178

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

FEATURES
This belt is a light and extremely comfortable back brace.
• Adjustable straps will help you find the ideal fit for you.
• Air-mesh fabric releases excess heat and moisture for all-day comfort.
• 4 Plastic stays provide firm grip, relief strain of lower back and protect weak muscles.
• Easy to care, dries quickly.

PERFORMANCE
Helps provide support to stiff, strained or aching back.

INDICATIONS
Lumbar Support for the orthotic fitting of the lower back.Recommended for the treatment of painful conditions of the lumbar spine. Lumbar sciatica, spondylosis, mild spondylosis, spondyloarthrosis, lumbar facet syndrome, lower back pain.

WARNINGS. PRECAUTIONS OF USE
Do not use the product if you are allergic to one of its components.
In case of any irritation occurrence, discontinue the use of the product and seek medical advice

SET UP AND OPERATION
Wrap the belt around your abdomen, make sure it is positioned on the centre of your back and secure the Velcro fasteners. Use the additional tension bands to adjust to the tightness that suits your own needs.

In addition to measuring the waist, the belt should be tried on to ensure a good fit and to understand how toput it on properly. Tighten the belt based on the intensity of the pain andtype of activity.

CARE INSTRUCTION
Hand wash in warm soapy water (maximum 30°C). Dry flat, away from sources of excessive heat (radiator, tumble dryer)

MEASUREMENT AND SIZE CHART (in Cm.)
The size is to choose in accordance with the waist circumference

S	M	L	XL	XXL	XXL
75 - 85	85 - 95	95 - 105	105 - 115	115 - 125	125 - 140

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