A-5071

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EN Wrist Support A-5071

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

FEATURES
This bandage will help support your wrist, reduce strain and discomfort, whilst allowing for normal daily functions.
Supports the wrist joint and helps recover from injury.

PERFORMANCE
Helps relieve tired wrist joints and muscles from swelling during activity.

INDICATIONS
For sequelae of fracture and sprains, heavy work, using vibrating instruments, sport activities. Stabilization of the wrist by adjustable closure

WARNINGS. PRECAUTIONS OF USE
Do not use the product if you are allergic to one of its components.
In case of any irritation occurrence, discontinue the use of the product and seek medical advice

SET UP AND OPERATION
Fit the bandage around your wrist and secure with the hook and loopfas ener.

CARE INSTRUCTION
Hand wash in warm soapy water (maximum 30°C). Dry flat, away from sour es of excessive heat (radiator, tumble dryer)

MEASUREMENT AND SIZE CHART (in Cm.)
The size is to choose in accordance with the wrist circumference

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